



U.S. DEPARTMENT OF
ENERGY

PNNL-27615

Radiation Detection and Nuclear Security Summer School

Activity 3A: Pulse Shape Discrimination and Fast Neutron Detection

J Baciak

June 2018



Pacific Northwest
NATIONAL LABORATORY

*Proudly Operated by **Battelle** Since 1965*

DISCLAIMER

This report was prepared as an account of work sponsored by an agency of the United States Government. Neither the United States Government nor any agency thereof, nor Battelle Memorial Institute, nor any of their employees, makes **any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights.** Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or any agency thereof, or Battelle Memorial Institute. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or any agency thereof.

PACIFIC NORTHWEST NATIONAL LABORATORY
operated by
BATTELLE
for the
UNITED STATES DEPARTMENT OF ENERGY
under Contract DE-AC05-76RL01830

Printed in the United States of America

Available to DOE and DOE contractors from
the Office of Scientific and Technical Information,
P.O. Box 62, Oak Ridge, TN 37831-0062

www.osti.gov
ph: (865) 576-8401
fox: (865) 576-5728
email: reports@osti.gov

Available to the public from the National Technical Information Service
5301 Shawnee Rd., Alexandria, VA 22312
ph: (800) 553-NTIS (6847)
or (703) 605-6000
email: info@ntis.gov
Online ordering: <http://www.ntis.gov>

Radiation Detection and Nuclear Security Summer School

Activity 3A: Pulse Shape Discrimination and Fast Neutron Detection

J Baciak

June 2018

Pacific Northwest National Laboratory
Richland, Washington 99352

1 Pulse shape discrimination and fast neutron detection

ADVANCE PREPARATION:

1. Read Knoll, pp. 569-582 (ch 16. fast neutron detection), 225-235 (ch 8. scintillation detectors). 700-701 (additional PSD)

THE NEUTRON SOURCE:

The source of neutrons used will be a ^{252}Cf source that contains about 1/3 of a microgram. 97% of decays are alpha decays while about 3% are spontaneous fission events emitted with an approximate Maxwellian/Watt spectrum. These neutrons are formed with a spread of kinetic energies ranging up to about 10 MeV (see lecture notes).

DETECTORS:

The detector we will use is a EJ-309 liquid scintillation detector. This contains a hydrocarbon liquid with a PMT, which is powered with a standalone negative high voltage supply. The digitizer is controlled by a computer program. The program is controlled with several options in an input script. You will have at your disposal a matlab script that can read in pulses and store them on your computer. That same script can also be used to calibrate the detector.

HIGH VOLTAGE: The scintillator PMTs should not exceed -2000 V.





Power supply (example) and Struck SIS3316-DT 16 Channel digitizer.

PROCEDURE:

EJ 309 Scintillation detector

(a) Calibration

We want to calibrate the detector to have a Compton edge of cesium 137 at 0.35 volts. This achieved using the matlab calibration script “calibration_read.m”. The same script will be used to read the various collected data sets.

The detector will be powered by negative high voltage from a standalone HV-power supply that is operated by touchscreen for up to 4 detectors (one will be used in this lab).

The detector output signal (“anode”) will be connected via BNC cable to a LEMO connector and then to the grey digitizer (any channel).

Using the “sis3316_test_new.exe” program, and the accompanying “sis3316.i” input script we shall collect a few pulses to verify that the detector is operating. At this stage we use a nearby Cs-137 source with the detector.

Once you see pulses take a data set of at least 20 000 pulses and run the “calibration_read.m” script to determine the current Compton edge value. Keep adjusting the voltage until a 0.35 V Compton edge is achieved (higher voltage will move it up and vice versa). What is a reasonably detection threshold in volts? What gamma-ray energy would that correspond to?

(b) Background measurement.

Take a background measurement of at least 10 000 pulses.

Where should we place the sources during this measurement? Why is the background rate so high?

(c) Measurements of neutron and gamma-ray radiation

1. Using an instructor supplied Cf-252 source take on the order of 100 000 pulses of combined neutron and gamma-ray data.
2. Using various gamma-ray energy sources take a measurement of 20 000+ gamma-ray pulses. (where/how should the neutron source be located/shielded at this time?)

3. Using shielding try to collect *ONLY* neutron data. What shielding should be used? At what location should the shielding be placed? (Collect 20 000+ pulses).

(d) instructor demonstration of tail-to-total pulse shape discrimination (PSD)

ANALYSIS:

1. Using the matlab script obtain the raw pulse waveforms from the 4 measurements (BG, neut+ γ , neut, γ).
2. Correcting for highly energetic cosmic rays in the data might be a good idea. Those pulses would be larger than the maximum recordable value of the digitizer and would then be distorted.
3. Apply a PSD algorithm to the waveforms, similarly to the instructor shown tail-to-total demonstration (has to be an algorithm other than the traditional tail-total integration comparison). You should see some separation in the combined data and can verify the PSD algorithm with the “pure” separate gamma-ray and neutron data sets. (Do you expect perfect agreement in either case)?
4. Explain your algorithm and its origin (created by you, obtained from references etc?)
5. Once you have a functioning PSD algorithm I want you to obtain the *background corrected* neutron and gamma-ray rate from the Cf-252 source. Note that there is a timing vector generated by the matlab script (“time_t”) which lists times of the each pulse in units of 4ns. The values are generated by the computer and will not restart with each measurement.
6. Estimate the fractional contribution to NSC background dose from gamma-rays vs fast neutrons assume no other radiations exist). You can assume a W_R -factor (used to be Q-factor) of 10 for neutrons and a W_R -factor of 1 for gamma-rays.



902 Battelle Boulevard
P.O. Box 999
Richland, WA 99352
1-888-375-PNNL (7665)

www.pnl.gov



U.S. DEPARTMENT OF
ENERGY