

PNNL-32645

Watchmen 3.1.0 Users Guide

Revision 7

February 2022

DT Keller
RE Wilson
JB Chapman
CJ McCall
CT Smith

DISCLAIMER

This report was prepared as an account of work sponsored by an agency of the United States Government. Neither the United States Government nor any agency thereof, nor Battelle Memorial Institute, nor any of their employees, makes **any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights.** Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or any agency thereof, or Battelle Memorial Institute. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or any agency thereof.

PACIFIC NORTHWEST NATIONAL LABORATORY
operated by
BATTELLE
for the
UNITED STATES DEPARTMENT OF ENERGY
under Contract DE-AC05-76RL01830

Printed in the United States of America

Available to DOE and DOE contractors from
the Office of Scientific and Technical Information,
P.O. Box 62, Oak Ridge, TN 37831-0062

www.osti.gov
ph: (865) 576-8401
fox: (865) 576-5728
email: reports@osti.gov

Available to the public from the National Technical Information Service
5301 Shawnee Rd., Alexandria, VA 22312
ph: (800) 553-NTIS (6847)
or (703) 605-6000
email: info@ntis.gov
Online ordering: <http://www.ntis.gov>

Watchmen 3.1.0 User Guide

Revision 7

DT Keller
RE Wilson
JB Chapman
CJ McCall
CT Smith

February 2022

Prepared for
the U.S. Department of Energy
under Contract DEAC0576RL01830

Pacific Northwest National Laboratory
Richland, Washington 99352

Revision Log

Watchmen 2 User Guide

Rev. No.	Date	Describe Changes	Pages Changed
0	5/06/2016	Original Issued	NA
1	7/22/2016	General Improvements	Throughout
2	1/11/2017	Additions: Status Indicators and notifications, RRR/ARR, Trending, Graphing, Calculations	Throughout
3	3/24/2017	Watchmen software fixes and enhancements, XECON calculations	Throughout
4	8/18/2017	Interactive Particulate	Throughout
5	9/16/2020	User authentication and user interface revamped	Throughout
6	9/15/2021	Performance Monitoring, Meta-Stable Interference analysis and multiple detector background analysis results	Throughout
7	1/25/2022	Additions: World, Calibration Dropdown Updates: Ad Hoc Query Tool, Performance Monitoring, Sample Review Panels, Stations	Throughout

Acronyms and Abbreviations

2-D	Two-dimensional
ARR	Automatic Radionuclide Report
CSV	comma-separated value
FWHM	full width at half maximum
ID	identification
IDC	International Data Centre
IMS	International Monitoring System
MDA	minimum detectable activity
MDC	minimum detectable concentration
MITER	Meta-stable Interference
PDF	Portable Document Format
PHD	pulse height data
PNNL	Pacific Northwest National Laboratory
QC	quality control
RN	radionuclide
ROI	region-of-interest
RRR	Reviewed Radionuclide Report
SAUNA	Swedish Automatic Unit for Noble gas Acquisition
SOH	state of health
scm	standard cubic centimeters

Contents

Revision Log.....	iii
Acronyms and Abbreviations	v
1.0 Watchmen.....	1.1
1.1 Product Functions.....	1.1
2.0 Viewer	2.2
2.1 Logging In.....	2.2
2.2 Dashboard	2.3
2.3 Menu	2.3
2.4 Review Queue	2.4
2.5 Thumbnails.....	2.5
3.0 Searching	3.1
3.1 Sample Review Panels	3.3
3.1.1 Sample Review at a Glance.....	3.3
3.1.2 XECON Report	3.6
3.1.3 Recalculate	3.7
3.1.4 Related Measurements	3.7
3.1.5 Calibration Dropdown.....	3.8
3.1.6 Automatic Radionuclide Report and Reviewed Radionuclide Report Views (currently not available)	3.8
3.1.7 Measurement Navigation Control	3.10
3.1.8 Status Indicators	3.10
3.1.9 Chart Features	3.12
3.1.10 Station Frequencies and Trends	3.13
3.1.11 All Spectra.....	3.14
3.1.12 Calibration Panel.....	3.15
3.1.13 Review Panel.....	3.16
3.2 Completing a Review – Alternative One	3.17
3.3 Completing a Review – Alternative Two.....	3.17
3.4 Selecting Isotope Checkboxes.....	3.17
3.5 Submitting the Review	3.18
3.6 QC Control Survey.....	3.18
4.0 Performance Monitoring.....	4.20
4.1 Performance Trends	4.20
4.2 Count Rate Control Trends	4.21
4.3 Gain Control Trends.....	4.22
4.4 Editing Control Limits	4.23
4.5 Creating Reports.....	4.24

5.0	Isotope Analyzer.....	5.1
6.0	Station Health	6.1
6.1	Station Health Application	6.1
7.0	Live View	7.1
8.0	SOH Overview	8.1
8.1	Station Performance Overview	8.1
8.2	Station Health.....	8.1
9.0	Ad Hoc Query Tool.....	9.1
10.0	World.....	10.4
11.0	Administration.....	11.5
11.1	Pending Reviews.....	11.5
11.2	Stations.....	11.6
11.3	Detector Configuration.....	11.12
11.4	Indicator Configuration.....	11.13
11.5	Users.....	11.1
12.0	Information.....	12.1
12.1	Server Information	12.2
12.2	Calculation Information	12.2
12.3	Legal.....	12.2
12.4	Notice	12.2
12.5	Developer Information	12.2

Figures

Figure 2-1. Hover to view sample details	2.4
Figure 3-1. Search Criteria Query	3.1
Figure 3-2. URL with Search Criteria Query	3.1
Figure 3-3. Search Results, Grid View	3.2
Figure 3-4. Search Results, List View	3.2
Figure 3-5. Search Results, Sample Review Panel	3.3
Figure 3-6. SAUNA Review	3.4
Figure 3-7. 2-D Histogram	3.5
Figure 3-8. 2-D Histogram Region Selection	3.6
Figure 3-9. XECON Report Button	3.6
Figure 3-10. XECON Report	3.7
Figure 3-11. Related Measurements	3.8
Figure 3-12. Calibration Dropdown with Custom Calibration	3.8
Figure 3-13. ARR/RRR details with MDC comparison outlined in blue	3.9
Figure 3-14. Measurement Navigation Control	3.10
Figure 3-15. Status Indicators	3.12
Figure 3-16. Measurement Information	3.12
Figure 3-17. Chart Features	3.13
Figure 3-18. Trend Panel	3.14
Figure 3-19. All Spectra Panel Beta-gamma Measurements (Normal/History View)	3.15
Figure 3-20. Calibration Panel	3.16
Figure 3-21. Review Panel	3.17
Figure 3-22. QC Control Survey	3.18
Figure 4-1. Performance Monitoring	4.20
Figure 4-2. Performance Trends tab	4.21
Figure 4-3. Count Rate Control Trends tab	4.22
Figure 4-4. Gain Control Trends tab	4.23
Figure 4-5. Editing control limits	4.24
Figure 4-6. Performance Monitoring Report	4.25
Figure 5-1. Isotope Analyzer	5.1
Figure 5-2. Isotope selector	5.1
Figure 5-3. Series parameters	5.2
Figure 5-4. Model selection and parameters	5.2
Figure 5-5. Isotope Analyzer plotting	5.3
Figure 5-6. Multi-isotope plot	5.3
Figure 6-1. Station Health Application	6.1

Figure 7-1. Live View Application	7.1
Figure 8-1. Station Performance	8.1
Figure 8-2. Station Health	8.2
Figure 8-3. Station Health Hover	8.2
Figure 9-1. Initial Ad Hoc Query Tool Page	9.1
Figure 9-2. Plotting a Single Variable	9.1
Figure 9-3. Plotting Multiple Variables	9.2
Figure 9-4. Adjusting the Date Range	9.2
Figure 9-5. Line Series Toggled	9.3
Figure 9-6. Moving Range Chart	9.3
Figure 10-1. World Page	10.4
Figure 11-1. Pending Sample Reviews Application	11.5
Figure 11-2. Station List	11.6
Figure 11-3. New Station Dialog	11.6
Figure 11-4. New Station Information Panel	11.7
Figure 11-5. Station Information Panel	11.8
Figure 11-6. Detector Calibrations	11.9
Figure 11-7. Viewing Calibrations	11.10
Figure 11-8. Station Trend Range Panel	11.11
Figure 11-9. Default Detector Configuration	11.12
Figure 11-10. Adding Custom Calibration	11.12
Figure 11-11. Mis-match Detector Dialog	11.12
Figure 11-12. Uploaded Detector Calibrations	11.13
Figure 11-13. Indicator Administration Start	11.13
Figure 11-14. Station Indicator Administration Configuration	11.14
Figure 11-15. Toggle E-mail Notifications	11.15
Figure 11-16. Indicator Configuration Detail	11.15
Figure 12-1. Information Detail	12.1

Tables

Table 2.1. Thumbnails	2.5
------------------------------------	-----

1.0 Watchmen

Watchmen is a research software application developed by Pacific Northwest National Laboratory (PNNL) that incorporates the scientific and operational expertise for reviewing data from treaty monitoring radionuclide stations. The radionuclide stations are part of the International Monitoring System (IMS), a worldwide network to monitor for nuclear explosions. The data the IMS produces are critical to determine if a radionuclide release event is from a nuclear explosion. Stations deliver their measurements and system status to the International Data Centre (IDC), which forwards it via email to all subscribers. Watchmen is capable of processing data from several radioxenon station types.

Screening of data in Watchmen may be done by several different users such as radionuclide analysts, evaluators, and data quality experts. This guide assists those users in navigating the application. The term Watchmen is used generically throughout this document to refer to any of the various components in the software application. The user interface that is viewed with a web browser is the primary focus of this user guide. Other components include a database to store measurements and state of health (SOH) data; and the data loader that monitors incoming emails, parses the data, populates the database, does the initial analysis, and routes data for review.

1.1 Product Functions

The Watchmen software is a toolset for monitoring station performance and diagnosing problems. It also provides basic screening of radioxenon measurements. In support of this goal, the software is equipped with tools for monitoring SOH data.

The primary functions of the software at the time of this report are:

- Graphical display of radionuclide measurements
- Data analysis to calculate estimated concentrations and minimum-detectable-concentrations (MDCs) for radionuclides
- Automated radioxenon
- Ability to trend measurement analysis results
- Ability to facilitate a review process for radionuclide measurements, including the ability to escalate
- Ability to override the station-provided energy gain calibration as part of data quality review processes for cases when the energy calibration is suspect
- SOH functionality, including ability to monitor specific status indicators of any station's health.
- Ability to monitor QC performance

2.0 Viewer

2.1 Logging In

Navigate in a web browser (Chrome or Firefox) to the Watchmen site. Ask an administrator for the URL. It will be in the form of `https://servername/webviewer2`, with *servername* as a placeholder for the name of the server or machine.

After authenticating as prompted (most likely by selecting a certificate from a smartcard) you will be prompted to request an account the first time you access the site.

Your account is either disabled or does not exist
[Request Access](#)


Click “Request Access” to create a disabled account in your name.

Your account is either disabled or does not exist
Notify the system administrator to activate your account, they will not be emailed.

Then, contact an administrator (Director of Operations) of Watchmen to have your account enabled.

They will not know you have requested access unless you inform them. Likewise, they will not be able to activate your account unless you click “Request Access”.


2.2 Dashboard



The first screen you will see is the dashboard which gives quick access to the most used tools in Watchmen.

2.3 Menu

Along the top of the screen is a menu.



This menu gives quick access to review your queue and all the tools available in Watchmen. If you are Director of Operations, you also see an Administrator option allowing configuration of stations and users.

2.4 Review Queue

You can access your own review queue and the queues for any roles you possess from the “Review Queue” menu or by selecting the user on the search page. The dashboard also includes a quick link to your queue.

The review queue holds all the measurements assigned to the user or role. Releasing or reassigning a measurement to another role removes it from the original queue.

Review Queue

Tools

My Queue

Lead Analyst

Senior Analyst

Director of Operations

Processing Engineer

USX01_002-2021/09/15-05:37

A S

USX01_002-2021/09/15-16:52

Q

16-17:37

A S

18-05:37

C S

Detector

USX01_002

Message Type

SAMPLEPHD

Reference ID

01202109141211G

Measurement ID

USX01_002-2021/09/15-05:37

Station

USX01

Collection Stop

2021-09-14T22:37:14.000Z

Acquisition Start






2021-09-15T05:37:15.000Z

Figure 2-1. Hover to view sample details

2.5 Thumbnails

On the review and search pages where the thumbnails are displayed, different types of files are color coded for a visual representation (Table 2.1).

Table 2.1. Thumbnails

Measurement Type	Color	Example Thumbnail
Sample	Purple	
Quality Control (QC)	Blue	
Gas Background	Green	
Detector Background	Yellow	
Spike/Calibration	Black	

3.0 Searching

Watchmen retrieves measurements from its own database; it is not possible to search for IMS measurements that have not been processed by Watchmen. Measurements can be searched by using any combination of the following criteria:

- Station
- Detector
- Type (measurement type)
- Acquisition Start FROM and TO date range
- Assigned to (for review)
- Quality (review quality)
- ID (Sample Reference Identification or Measurement ID)
- Preliminaries (toggle to set to exclude, only “FULL” measurements are returned)

NOTE: Quality and Assigned To are mutually exclusive because once reviewed and given a quality, it is no longer assigned to anyone.

The primary set of criteria are “Station,” “Detector,” “Type,” “Acquisition Start” and “Assigned To,” which are immediately visible upon arriving at the search page. The second set of criteria are “Quality,” “Category,” “ID,” and “Preliminaries,” which can be displayed by clicking ellipses button (⋮) at the upper right. When searching by “ID”, the operator may use the Sample Reference ID (e.g., 0188888882342G) or the Measurement ID (e.g., ABC01_001-2016/04/08-15:16). Searching on these IDs may have multiple results if the operator toggles “Preliminaries” to included (default is excluded) or if the IMS retransmitted an altered measurement.




Figure 3-1. Search Criteria Query

As criteria is entered, the results are filter. Once you have selected a measurement, the address bar of the browser will also include the Watchmen-generated ID that uniquely identifies the measurement in the database. Sending someone the URL in the address bar of the measurement you are viewing will allow them to see the same measurement and search criteria, without searching.




Figure 3-2. URL with Search Criteria Query




Figure 3-3. Search Results, Grid View

There is an optional list view (Figure 3.4) that is used for reassigning samples, that can be viewed by pressing the tab icon in the upper right corner, just below the ellipse button. A check box is located next to each sample for selection. To select multiple contiguous samples: 1) press and hold shift and then 2) click the first and last sample in the range. If a sample has already been reviewed, the check box will not appear and cannot be selected, including when a range of samples is selected.

Station

USK77

Detector

Search...

Type

Sample

Acquisition Start

Start date

End date

Assigned to

Lead Analyst

Sort By Date

Oldest

Newest

Quality

Search...

Category

Search...

ID

Search...

Preliminaries

Excluded

<input type="checkbox"/>	ID	Station ID	Measurement ID	Type	Collection Stop	Quality	Assigned To	Reassign
<input type="checkbox"/>	1794076001	USK77	USK77_001-2016/01/21-10-12	SAMPLEPHD	01/20/2016 19:14:00		Lead Analyst	
<input type="checkbox"/>	1796292001	USK77	USK77_002-2016/01/21-22-12	SAMPLEPHD	01/21/2016 07:14:00		Lead Analyst	
<input type="checkbox"/>	1798588001	USK77	USK77_001-2016/01/22-10-12	SAMPLEPHD	01/21/2016 19:14:00		Lead Analyst	
<input type="checkbox"/>	1800752001	USK77	USK77_002-2016/01/22-22-12	SAMPLEPHD	01/22/2016 07:14:00		Lead Analyst	
<input type="checkbox"/>	1802868001	USK77	USK77_001-2016/01/23-10-12	SAMPLEPHD	01/22/2016 19:14:00		Lead Analyst	
<input type="checkbox"/>	1805036001	USK77	USK77_002-2016/01/23-22-12	SAMPLEPHD	01/23/2016 07:14:00		Lead Analyst	
<input type="checkbox"/>	1807208001	USK77	USK77_001-2016/01/24-10-12	SAMPLEPHD	01/23/2016 19:14:00		Lead Analyst	
<input type="checkbox"/>	1809456001	USK77	USK77_002-2016/01/24-22-12	SAMPLEPHD	01/24/2016 07:14:00		Lead Analyst	
<input type="checkbox"/>	1811648001	USK77	USK77_001-2016/01/25-10-12	SAMPLEPHD	01/24/2016 19:14:00		Lead Analyst	
<input type="checkbox"/>	1813864001	USK77	USK77_002-2016/01/25-22-12	SAMPLEPHD	01/25/2016 07:14:00		Lead Analyst	
<input type="checkbox"/>	1816108001	USK77	USK77_001-2016/01/26-10-12	SAMPLEPHD	01/25/2016 19:14:00		Lead Analyst	
<input type="checkbox"/>	1819428001	USK77	USK77_002-2016/01/26-22-12	SAMPLEPHD	01/26/2016 07:14:00		Lead Analyst	
<input type="checkbox"/>	1820492001	USK77	USK77_001-2016/01/27-10-12	SAMPLEPHD	01/26/2016 19:14:00		Lead Analyst	
<input type="checkbox"/>	1822644001	USK77	USK77_002-2016/01/27-22-12	SAMPLEPHD	01/27/2016 07:14:00		Lead Analyst	
<input type="checkbox"/>	1824824001	USK77	USK77_001-2016/01/28-10-12	SAMPLEPHD	01/27/2016 19:14:00		Lead Analyst	
<input type="checkbox"/>	1827064001	USK77	USK77_002-2016/01/28-22-12	SAMPLEPHD	01/28/2016 07:14:00		Lead Analyst	
<input type="checkbox"/>	1829964001	USK77	USK77_001-2016/01/29-10-12	SAMPLEPHD	01/28/2016 19:14:00		Lead Analyst	

Figure 3-4. Search Results, List View

Once the desired measurement has been found, click on the thumbnail to view the Sample Review Panels (see Section 3.1 for details).




Figure 3-5. Search Results, Sample Review Panel

3.1 Sample Review Panels

3.1.1 Sample Review at a Glance

Selecting a sample with a single left-click of the mouse in the search will bring up the Sample Review Panel (Figure 3.6). The interface will display the Sample Review Panel, for the selected measurement, below the graphical list of searched measurements. This panel provides detailed information (e.g., activity concentration, measurement times, etc.) on the currently selected sample along with the associated detector and gas background measurements, if the sample is from a beta-gamma type station. The measurement review layout displays the Detector name and Acquisition Start date at the top along with action buttons. The action buttons include Xeon Report (Section 3.1.2), Delete, Recalculate (Section 3.1.3); Related Measurements (Section 3.1.4); and Measurement Navigation Control (Section 3.1.6).

WARNING: The Delete button will delete the measurement and its analysis from the database.¹

Next, the Overview panel contains the status indicators, Spectrum chart, two-dimensional (2-D) histogram (if a beta-gamma sample), a table with information about the measurement, and log entries.

¹ Data can only be recovered by locating the original file in the archive and placing the file in the incoming directory for re-ingesting and re-analysis by Mothman.

Following the Overview is the Trend panel (if a beta-gamma sample), the All Spectra panel (if a beta-gamma sample), Calibration panel, and finally the Review panel.

While the interface contains a similar layout for all measurements, the populated information will vary depending on the station and calculation type. Figure 3.6 and Figure 3.7 show a Swedish Automatic Unit for Noble Gas Acquisition (SAUNA) sample for the screen layout.




Figure 3-6. SAUNA Review

The overview contains the 2-D histogram with beta energy (in keV) on the x-axis and the gamma energy (in keV) on the y-axis (Figure 3.7). Regions of interest (ROIs) are outlined with the rectangular bars on the graph. The ROIs can be turned off or on by clicking the 'Hide Regions' and 'Show Regions' button respectively. The box is bolded when regions are turned on.




Figure 3-7. **2-D Histogram**

Individual ROIs can be turned off or on by clicking the filter icon (which is only visible when the 'Hide Regions' option is available) and selection them from the list regions. If only some of the ROIs are being displayed, and the user clicks 'Hide Regions', the subsequent use of the 'Show Regions' button will restore all the regions (does not preserve previous filtering).




Figure 3-8. 2-D Histogram Region Selection

The ROIs are used to compute the four reported radioxenon isotopes ($^{131\text{m}}\text{Xe}$, $^{133\text{m}}\text{Xe}$, ^{133}Xe , and ^{135}Xe).

3.1.2 XECON Report

You can perform XECON calculations and generate a textual report on your selected sample. Click on the “Xecon Report” button (Figure 3.8) and a new browser tab will appear. This new tab will produce the XECON Report (Figure 3.9).

NOTE: XECON version 1.2.1 is used by the software.




Figure 3-9. XECON Report Button

Xecon Report for USX79_004, 2022/01/21 (SAMPLEPHD)

XENON SAMPLE ANALYSIS			
Created by:xecon 1.2.1 (2010-05-17)			
Sample:			
/analysis_tools/xecon/tmp/sample.phd			
Gas background:			
/analysis_tools/xecon/tmp/gasbk.phd			
Detector background:			
/analysis_tools/xecon/tmp/detbk.phd			
Collection start: 2021/11/25 12:02:16.0 Collection stop: 2021/11/25 18:02:16.0			
Air volume: 29.34 m3			
Collection time: 21600 sec			
Processing time: 9756 sec			
Xenon volume: 2.552230 +/- 0.038280 cm3			
Sample acq start: 2021/11/25 20:44:52.0			
Sample acq times: 43141.300000 (real) 43260.000000 (live)			
Gasbk acq start: 2021/11/25 08:54:36.0			
Gasbk acq times: 38021.000000 (real) 38085.000000 (live)			
Detbg acq start: 2021/09/09 16:48:32.0			
Detbg acq times: 254929.000000 (real) 259262.000000 (live)			
Det subf sample: 0.169229			
Det subf gasbk: 0.149143			
Time difference between s and g measurements: 42616 sec			
Gas background subtraction factors:			
F_133: 1.058912 F_135: 0.441220 F_131m: 1.100706 F_133m: 0.961737			
Det subf gasbk: 0.149143			
Energy calibration tweaking parameters (offset gamma, tweak gamma, offset beta, tweak beta): 0.000000, 1.000000, 0.000000, 1.000000			
Memory effect correction: 1.000000			
Gross counts in ROI:s:			
ROI	Sample	Gasbk	Detbk

1	165	146	776

Figure 3-10. XECON Report

3.1.3 Recalculate

If a detector's energy coefficients have been changed (see Stations Section 9.2), the new calculations are not automatically applied to measurements that were already in your review queue, as new coefficients are only applied to samples that are inserted to the database after the coefficients were applied, or if you click Recalculate. If the coefficients have been changed recently, select the Recalculate button located on the right, in line with the Overview panel to apply the new calculations. Refresh the screen when the calculations have finished to view the modified data.

NOTE: This change is permanent, and the measurement will be stored in the database with the revised calculations. In addition, when the calculations are rerun, the sample will be assigned back to the default reviewer for the station, which may be a different review queue if the measurement had been assigned to a lead analyst for review.

3.1.4 Related Measurements

Once you have selected a measurement, you can navigate to any of its related measurements (Backgrounds, Quality Checks, or Preliminaries) by using the Related Measurements drop-down menu located on the right side of the screen (see Figure 3.10). If there are no related measurements in the

database, the drop-down menu will be disabled. Selecting a measurement in the Related Samples list will take you to the selected measurement.




Figure 3-11. Related Measurements

3.1.5 Calibration Dropdown

The calibration dropdown can be used to view the sample with a calculation from a custom calibration. The ROIs, Indicators, Isotopes, Trend Panel, and Calibration Panel all update when switching to a custom calibration. The calibration is disabled in the dropdown if there is no calculation for the calibration. If a calibration is disabled, the ‘Recalculate’ button can be used to generate a calculation for the calibration.



Figure 3-12. Calibration Dropdown with Custom Calibration

3.1.6 Automatic Radionuclide Report and Reviewed Radionuclide Report Views (currently not available)

Automatic Radionuclide Report (ARR) and Reviewed Radionuclide Report (RRR) views are available from the Related Measurements menu. If an ARR or RRR is available, it will be clickable. Otherwise, they will be grayed out if they do not exist. Once you click on the ARR or RRR (see Figure 3.10), the web browser window will navigate to the detailed ARR or RRR data.

The detailed ARR or RRR displays the report information along with a percent difference comparison to the Watchmen MDC calculations. If the calculation difference is greater than 5%, the difference is highlighted in red. You can also export the MDC calculations to comma-separated value (CSV) format—indicated by the “Download Summary” link. See Figure 3.11—the MDC comparison section is outlined in blue.

Radionuclide Report Comparison - ARR

Creation Date: 2016/12/26 17:51:09 UTC
Sample Arrival Time: 2016/12/26 17:40:03 UTC
Time difference from receipt of raw data to report creation: 0d 0h 11m 6s

Sample Information

Station ID:	USX74	Detector Code:	USX74_006
Authenticated:	NO		
Station Location:	Ashland, KS, USA		
Detector Description:	Detector #6 for USX74 Noble Gas Test Station.		
System Technology:	SAUNA		
Sample Reference ID:	74201612251111G		
Sample ID:	3907557		
Stable Xe Volume:	1.07 mL	Sample Type:	Gas
Collection Start:	2016/12/25 11:19:41 UTC	Sampling Time:	0d 12h 0m 2s
Collection Stop:	2016/12/25 23:19:43 UTC	Processing Time:	0d 7h 13m 42s
Acquisition Start:	2016/12/26 06:33:25 UTC	Acquisition Time:	0d 10h 58m 17s
Acquisition Stop:	2016/12/26 17:31:42 UTC		

Measurement Categorization

Categorization Legend

Level A	=	Clean spectrum - No Xenon is present in the sample.
Level B	=	Xenon detection within the typical range for the station.
Level C	=	Anomalous Xenon detection.

Isotope Category

Isotope	Nuclide detected	Abnormal_limit (mBq/m3)	Category
131m Xe	NO	1.61e-1	A
133m Xe	NO	1.24e-1	A
133 Xe	NO	3.70e-1	A
135 Xe	NO	9.62e-1	A

Spectrum Category: A

Activity Summary and Minimum Detectable Concentration for Xenon Isotopes

Radon Level in Xenon sample

Nuclide	Half-Life	Area	%RelErr
---------	-----------	------	---------

Radon counts in Xenon sample: 93

Xenon Isotopes

Net Count Calculation Analysis Method

Beowulf Tolerance Threshold: 5%

Download summary: [\[X\]](#)

Nuclide	Half-Life	LC	Conc (mBq/m3)			MDC			%RelErr
			IDC	Beowulf	% Diff	IDC	Beowulf	% Diff	
131m Xe	1.03e+6	4.00e-2	< LC	5.13e-2	N/A	1.00e-1	1.54e-1	5.36e+1	N/A
133m Xe	1.89e+5	4.00e-2	< LC	5.30e-2	N/A	1.10e-1	1.30e-1	1.85e+1	N/A
133 Xe	4.53e+5	9.00e-2	< LC	8.77e-3	N/A	1.90e-1	2.07e-1	9.07	N/A
135 Xe	3.29e+4	3.90e-1	< LC	9.52e-1	N/A	8.20e-1	7.33e-1	1.06e+1	N/A

Processing Specific Parameters and Results

Method 1 (Net Count Calculation)

ROI Net Count Results

ROI	Nuclide	Net Counts	Abs Net Error	LC	Efficiency	Abs Eff Error
1	214 Pb	2.17e+1	1.04e+2	1.49e+1	N/A	N/A

Figure 3-13. ARR/RRR details with MDC comparison outlined in blue

3.1.7 Measurement Navigation Control

Once you have selected a measurement, it is possible to navigate to measurements on the same station or detector using the Measurement Navigation Control located on the right side of the screen (see Figure 3.12). You can specify if you would like to navigate based on the Station or Detector (the leftmost drop-down menu) as well as the type of measurement (rightmost drop-down menu). Arrows on either side of the navigation control will navigate to the previous (leftmost arrow) or next (rightmost arrow) measurement chronologically when clicked. If there is no measurement available, the arrows are disabled.




Figure 3-14. Measurement Navigation Control

3.1.8 Status Indicators

Multiple status indicators, such as xenon volume (ccXe), acquisition time and data arrival time, xenon activity and isotope ratios, and the presence of related measurements, are displayed on the top of the Overview section in the Sample Review Panel. The indicators are categorized by Operation, Values, and Data Integrity. These status indicators display as red, yellow, or green. If they are green, they are collapsed by default. To expand collapsed indicators, click on the category. Hover over the indicator with your mouse to display a pop-up with additional details such as bounds for the parameter. The concentrations listed in the overview panel also contain color-coded information that corresponds to the status indicators.

- Data Integrity
 - Analysis
 - Red – Calculations are missing.
 - Green – Calculations are completed.
 - Detector Background
 - Red – A detector background file is not present.
 - Yellow – A detector background file is outdated.
 - Green – A detector background file is present.
- IDC
 - Gas Background (on compatible station types)
 - Red – A gas background file is not present.

- Yellow – A gas background is outdated.
 - Green – A gas background file is present.
 - Quality Control – The quality control (QC) status lights for the Sample files are controlled by reviewing the corresponding QC file (Figure 3.13).
 - Green – The corresponding QC files were present and reviewed as good.
 - Yellow – The corresponding QC files were present and have not been reviewed yet or were not marked as a good file.
 - Red – All other conditions, including the corresponding QC files are not present.
 - Report Time
 - Data Arrival - Within the time specified by the indicator
- Operation
 - ccXe – Cubic centimeters of xenon (setpoint is based on SAUNA ideal of 1 sccm/24hr)
 - Red – Volume of xenon is greater than 0.5 sccm from the set point
 - Yellow – Volume of xenon is between 0.25-05 sccm from the setpoint
 - Green – Volume of xenon is within 0.25 sccm of setpoint
 - Collection
 - Red – Difference from expected collection time is greater than 2%
 - Yellow – Difference from expected collection time is between 1-2%
 - Green – Difference from expected collection time is less than 1%
 - Processing
 - Red – Difference from expected processing time is greater than 2%
 - Yellow – Difference from expected processing time is between 1-2%
 - Green – Difference from expected processing time is less than 1%
 - Acquisition Real Time - Within the tolerance for the station and sample type
- Values
 - Xe133 – Red, yellow, green to show Above, Above Half, or Below MDC (respectively)
 - Xe135 – Red, yellow, green to show Above, Above Half, or Below MDC (respectively)
 - Xe131m – Red, yellow, green to show Above, Above Half, or Below MDC (respectively)
 - Xe133m – Red, yellow, green to show Above, Above Half, or Below MDC (respectively)
 - 131m/133 Ratio – Red, yellow, green to show if the ratio of the isotopes is below the minimum detectable concentration. Additional information about concentration displayed.
 - 135/133 Ratio – Red, green to show if the ratio of the isotopes is below the minimum detectable concentration. Additional information about concentration is displayed.

- Gas – Detector – Detects levels of counts and coincidence.
- Measurement information (Figure 3.14) – The metadata associated with the file collection such as acquisition time and station. The IDs for the corresponding gas and detector background files are also displayed. This information is displayed in the scrollable pane.




Figure 3-15. Status Indicators




Figure 3-16. Measurement Information

3.1.9 Chart Features

The charts throughout the interface have a built-in zoom feature. With the mouse cursor over an area of interest in a chart, hold the ALT key and mouse-scroll to zoom in. You can hold ALT+SHIFT and mouse-scroll to zoom only the x-axis, or ALT+CTRL and mouse-scroll to zoom only in the y-axis. You also can click-and-drag to zoom into an area of interest. Once zoomed in, you can right-click to zoom

back out. Points in trend graphs (Figure 3.15) are clickable and linked to the measurement. Clicking the point will navigate to the sample measured.

The small download icon (📄) may be used to download the data to a CSV format. Select the download icon, then use the drop-down menu for Export to CSV.

Click on the y-axis to toggle between log and linear.

Each chart legend controls which data are displayed on the chart. Click on an item in the chart legend to turn the display of that specific data on and off.




Figure 3-17. Chart Features

3.1.9.1 Chart Panels

There are two detailed charting panels, Trend, and All Spectra. After selecting a measurement, the data from the different charts may be viewed by selecting the different tabs within these panels.

3.1.10 Station Frequencies and Trends

Two different displays for station frequencies and trends are supported in the Watchmen viewer - the frequency histogram and the trend over time plot as seen in Figure 3.16. The frequency histogram shows the frequency of the activity concentrations for each of the four radioxenon isotopes as well as the volumes of xenon for the samples. The trend data provide a different presentation of the same data to illustrate trends in the data over time. Both the frequency and trend data are shown together.

Two methods exist to display the frequency and trend data from a station. Firstly, the trends data is displayed along with the measurement information automatically (if available) on the search and review screens. The second is in the Station application and is described in Section 9.2.




Figure 3-18. Trend Panel

3.1.11 All Spectra

The All Spectra panel (Figure 3.17) shows different views of the selected measurement, indicated by the tabs. You can toggle the Normal and History tabs. Normal shows current measurement data from all relevant measurements. History shows aggregated data from the detector over one day, one week, and one month. See Section 3.1.10 for additional information on the chart features. The All Spectra panel contains the following tabs of graphs:

- Gamma Coincidence Spectrum: only gamma energy values where the gamma and beta spectra are in coincidence are displayed on the coincidence beta spectrum.
- Beta Coincidence Spectrum: only beta energy values where the beta and gamma spectra are in coincidence are displayed on the coincidence beta spectrum.
- Region 4 – Beta Coincidence: 30-keV Gamma, Coincidence Beta Spectrum: the plot is similar to the coincidence beta spectrum chart previously described; however, only coincidence data in the 30keV gamma region of interest are displayed.
- Gamma Singles Spectrum: The entire detector response for the gamma detector is displayed on the singles gamma spectrum. All values, including background are displayed without a coincidence being required.
- Beta Singles Spectrum: The entire detector response for the beta detector is displayed on the singles beta spectrum. All values, including background are displayed without a coincidence being required.




Figure 3-19. All Spectra Panel Beta-gamma Measurements (Normal/History View)

3.1.12 Calibration Panel

The values from the header for each file are displayed on the calibration panel (Figure 3.18). Select the type of values from the tabs.

- ROI Limits – The region of interest (ROI) limits for beta and gamma
- Ratios – The interference ratios for each region of interest
- Gamma
 - Efficiency - The gamma energy, efficiency, and efficiency uncertainty
 - Energies – The gamma energy, channel, and channel uncertainty
 - Resolutions – The gamma energy, full width at half maximum (FWHM), and FWHM uncertainty
- Beta
 - Energies - The beta energy, channel, channel uncertainty, and decay mode
 - Resolution - Beta energy, FWHM, and FWHM Uncertainty
- Beta Gamma Efficiency – ROI, nuclide name, and efficiency with uncertainty.

Calibration

ROI Limits Ratios Gamma Beta Beta Gamma Efficiency Software

ROI	Beta Start	Beta Stop	Gamma Start	Gamma Stop
1	10.7817	617.203	319.271	380.37
2	10.7817	921.752	214.135	279.124
3	10.7817	357.285	65.9624	98.3807
4	10.7817	390.203	20.509	41.8308
5	82.1945	151.89	20.509	41.8308
6	171.279	251.844	20.509	41.8308
7	10.7817	74.372	20.509	41.8308
8	260.786	390.203	20.509	41.8308
9	171.279	393.549	20.509	41.8308
10	10.7817	151.89	20.509	41.8308

Figure 3-20. Calibration Panel

3.1.13 Review Panel

Review and promote functions enable analysts to read and submit reviews and promote measurements for review by others. To complete the review on a sample, find the review panel at the bottom of the selected measurement user interface. You can also press <AccessKey>+R to jump to the review panel, <AccessKey> is typically Alt, Alt+Shift, or Ctrl+Alt, depending on the browser used. The review panel is the last panel in all measurement user interfaces. This panel contains the interface shown in Figure 3.19.

Use the “Send To” drop-down menu to assign the measurement to another user, or the “Quality” drop-down menu to complete a review.

The process of sample review is intended to be completed in chronological order. If it is completed in this order, then the QC file corresponding to the Sample file will always be reviewed first, thus changing the status lights to the expected status when reviewing the Sample file. The QC review should include a check of the energy calibration and acquisition time to verify that the results were as expected. See Section 3.1.14 of this guide for additional information on modifying the energy calibration if the QC review reveals results that need to be modified.

Send To

Quality

Good

Comment Templates

Comment

Enter text or choose from the templates

Isotopes to Report

ccXe

1.05 ± 0.10

	Isotope	Concentration	MDC
<input checked="" type="checkbox"/>	131m Xe	-0.021 ± 0.021 $\frac{\text{mBq}}{\text{m}^3}$	0.096 $\frac{\text{mBq}}{\text{m}^3}$
<input checked="" type="checkbox"/>	133 Xe	0.062 ± 0.037 $\frac{\text{mBq}}{\text{m}^3}$	0.13 $\frac{\text{mBq}}{\text{m}^3}$
<input type="checkbox"/>	133m Xe	0.0015 ± 0.024 $\frac{\text{mBq}}{\text{m}^3}$	0.097 $\frac{\text{mBq}}{\text{m}^3}$
<input type="checkbox"/>	135 Xe	0.25 ± 0.15 $\frac{\text{mBq}}{\text{m}^3}$	0.50 $\frac{\text{mBq}}{\text{m}^3}$

Submit

Report to be generated

US NATIONAL DATA CENTER GENERATED REPORT

REVIEWED RADIONUCLIDE REPORT

Xenon Version

=====

SAMPLE INFORMATION

=====

Station ID:

USX77

Detector ID:

USX77_004

Station Location:

Wake Island 19.3 N 166.6 E

Sample Quantity:

1.0496 cc Xe

Sample Type:

Gaseous

Collection Stop:

1/23/22, 1:01:47 PM GMT

This is a BetaGamma PINL BetaGamma 7 ROI analysis

=====

ACTIVITY SUMMARY

=====

FISSION-PRODUCT RADIOACTIVITY:

Isotope	MDC (mBq/SCM)	Concentration (mBq/SCM)	Error (mBq/SCM)
XE131m	0.0958	-0.0214	0.0214
XE133	0.1274	0.0615	0.0374

Figure 3-21. Review Panel

3.2 Completing a Review – Alternative One

One alternative is to pass the measurement to another user for review using the choices in the Send To drop-down menu. This drop-down menu contains the roles of “Lead Analyst,” “Senior Analyst,” “Process Engineer,” and “Director of Operations.” Users are mapped to the roles of “Lead Analyst,” “Senior Analyst,” etc., through the “Manage Users” application (see Section 9.3). After selecting the role from the drop-down menu, click the “Submit” button to transfer the spectra from your queue to the selected role.

3.3 Completing a Review – Alternative Two

The second alternative for a review is to complete the review by selecting a choice from the “Quality” drop-down menu of “Good,” “Warn,” or “Bad.”

A comment may be entered for either passing spectra or completing the quality review. The comments and results of the review are logged and subsequently displayed in the log. Some frequently used comments have been included in the “Comment Templates” drop-down menu. Select one of these if desired and then additional text may be entered.

3.4 Selecting Isotope Checkboxes

The Isotope checkboxes at the bottom of the Review page may be pre-selected depending on the calculations. These selections may be changed by checking or unchecking the box. If a checkbox is selected, the Reviewed Radionuclide Report will be generated for preview, and saved to the file system

upon submission. If no checkboxes are selected, the comment and action status will still be logged, but no Reviewed Radionuclide Report will be saved.

Checkboxes appear in the Isotopes to Report section once a Quality selection is picked. These checkboxes may or may not appear, depending on the measurement type.

3.5 Submitting the Review

After completing the review options, click the Submit button to complete the review.

3.6 QC Control Survey

The analysts reviewing the files can determine if they need to modify the energy calibration to change the gain in cases where the energy of peaks do not line up with known energy values (Figure 3.20).




Figure 3-22. QC Control Survey

Normally when the QC file is in the review queue and has been opened from the review queue, the analyst will verify that the Log Entries shows the automatic peak fitting was successful (log displays QC Peaks & Coefficients Calculated as shown in Figure 3.20), and that the peaks are as expected in the plot. The analyst will then click the Review Sample link to review the QC files in a process similar to the sample files.

The trend chart located at the bottom of the QC Review Normal page displays the peak centroids, area, and resolution over time for the QC files. A more comprehensive analysis may be accomplished by utilizing the Performance Monitoring tool.

4.0 Performance Monitoring

The Performance Monitoring module provides a tool for charting data collected from a single station-detector pair over long periods of time. To select data, analysts must select a detector, a file type, and a data range. Analysts may optionally select a station. If a station is selected, the list of detectors is limited to detectors associated with that station. When viewing Quality Check data, the cesium-137 decay correction feature is available and enabled by default. If correction is not desired, such as when viewing data from a system which uses isotopes other than cesium-137, it can be disabled using the toggle switch. The Get Data button will become enabled when all required fields are filled out. Clicking it will fetch the data with the given parameters and display it on a series of charts.




Figure 4-1. Performance Monitoring

After the data has loaded, it is displayed on several charts across three tabs. All the charts share several common features: Each chart supports zooming and panning as well as hiding of individual series (see section 3.1.8 for details). Clicking on a data point will open the measurement associated with that data point in a new tab. The data on each chart can be filtered. The Raw filter displays all selected data. The Seven Day Average plots the mean of each seven day period within the selected time range. And the Smooth Filter smooths the selected data using the Savitzky-Golay filtering algorithm.

4.1 Performance Trends

The Performance Trends tab is the default tab. It contains five charts: Coincidence, Singles, Compton Scatter Intercepts, Compton Scatter Intercept Resolution, and Channel Coefficients. If there is no data to display on a given chart, it will be inactive and display the message “No data found”. This feature can be seen on the Channel Coefficients chart in the screenshot below.




Figure 4-2. Performance Trends tab

4.2 Count Rate Control Trends

The Count Rate Control Trends tab displays one chart for beta total singles, one for gamma total singles, and then one coincidence chart for each ROI. Each chart displays a single series. While these charts display data that is also on the Performance Trends tab, they also control limits. The default control limits are ± 2 standard deviations (2σ) and ± 3 standard deviations (3σ). Lines at 2σ are orange, while those at 3σ are red. When data points violate one or more limits, a warning is generated, and the background of that data point is highlighted with the color of the highest priority limit violated. In the screen shot below, the data point violated the 2σ limit but not the 3σ limit, so its background is highlighted orange. The total number of warnings for all charts on this tab is indicated next to the tab's name.




Figure 4-3. Count Rate Control Trends tab

4.3 Gain Control Trends

The Gain Control Trends tab displays charts for beta and gamma Compton scatter intercepts and beta and gamma Compton scatter intercept resolutions. Each chart displays just one data series. The Compton scatter intercepts charts also have control limits. By default, these are static values. Like the Gain Control Trends tab, when data points violate one or more limits, a warning is generated, and the background of that data point is highlighted with the color of the highest priority limit violated. In the screen shot below, the data has violated the more extreme limit, so the warning is red. Note that when adjacent data points both exceed the same limit(s), only one warning is generated. Thus, the chart below contains only one warning. The total number of warnings for all charts on this tab is indicated next to the tab's name.





Figure 4-4. Gain Control Trends tab

4.4 Editing Control Limits

Details for the control limits on the Count Rate Control Trends and Gain Control Trends tabs can be viewed by clicking the Edit Control Limits button. The resulting table is populated with the control limits from the currently selected tab. Thus, to view control limits associated with charts on the Count Rate Control Trends tab, one must have the Count Rate Control Trends tab selected.

To create a new limit, select an existing control limit in the table or enter values into the control limits editor without selecting a control limit. When done, click Add to add the new limit to the table. To update a limit, select it in the table, modify any values, and click Update. An individual control limit can be deleted by clicking the  icon in its row. After changing the control limits, the charts and the warnings displayed on them will be updated. Note that any changes to control limits are *not* persistent. When the page is refreshed or new data is requested, any changes made to control limits are lost. If you would like to go back to the default limits, simply click the Get Data button.

In the screenshot below, a custom warning has been added to the beta total singles chart.




Figure 4-5. Editing control limits

4.5 Creating Reports

To create a monitoring report which contains the current data, click the Create Report button. This will open a new page with the same data in the same charts, but the layout is optimized for printing. Note that changes made to control limits *are* reflected on the monitoring report page. The parameters used to select the data are always included in the report. Additional notes can be entered in the Report Notes box using the Markdown language. When done reviewing and annotating, the report can be printed to a hard copy or a PDF using the browser's built-in print feature.

Report Notes

Notes can be formatted using the Markdown language.

[Learn more about Markdown](#)

Monitoring Report

Created on 2021-08-25T23:20:54.447Z

Notes

- Station: USX75
- Detector: USX75_008
- Type: QCPHD
- Start date: 2021-07-15T00:00:00.000Z
- End date: 2021-07-27T00:00:00.000Z




Figure 4-6. Performance Monitoring Report

5.0 Isotope Analyzer

The Isotope Analyzer module provides a charting tool for RN Evaluators to compare sample data with data from known sources or models, to determine or rule out possible sources of the radionuclide release. The Isotope Analyzer module allows the evaluators to select models and specify model parameters, select sample data of interest, and to configure one or more charts displaying model and sample data. A number of graphing options are available, including plotting isotopes and isotope ratios with respect to time, plotting ratios of isotopes (triple isotope plots), selecting the isotopes produced at a number of hold up times and irradiation times for different types of nuclear reactors and medical isotope production facilities, and allowing the user to add relevant samples collected on particulate and radionuclide systems. This provides a graphical comparison of real data sets to data sets from models.

The screenshot shows the 'Isotope Analyzer' interface. It has four main sections: 'Isotopes of Interest' with two empty dropdown menus; 'Sample Data' with a green '+ add series' button; 'Models' with a table of input fields (Model Type, Irradiation Time, Separation Time, Scalar, Start Date-Time) and a green '+ create chart' button below it.

Figure 5-1. Isotope Analyzer

The first step to using Isotope Analyzer is choosing the isotope(s) of interest. The chosen isotope(s) will be available for plotting on the charts and will have a column when entering data. Isotopes can be added or removed at any time. Removing an isotope will delete all hand entered activities for that isotope.

This screenshot shows the 'Isotope Analyzer' interface with the 'Isotopes of Interest' section active. A dropdown menu is open, showing a list of isotopes from 127 to 137. The isotope 135 is highlighted in blue. The background shows the same interface as Figure 5-1, but with the dropdown menu obscuring part of the 'Sample Data' and 'Models' sections.

Figure 5-2. Isotope selector

After choosing the isotopes you are interested in, data from sample measurements can be entered as a series into the Sample Data panel. Data from separate locations or equipment can be entered as different series. If you have multiple measurements from the same piece of equipment you can enter each and record the time of the measurement.

Values for measurements may be entered as activities or concentrations. The same unit must be used throughout all measurements or an incorrect conclusion will result. The “Dilution” of the entered values can be adjusted to normalize the entered values. The dilution value is a multiplicative term which will scale up or down the values. An example of the use of dilution would be to normalize activities measurements to concentration values by accounting for sample size.

Sample Data

Series Name:

Date-Time	Dilution	133 Xe	131m Xe	133m Xe
<input type="text" value="2015-07-01 00:00:00"/>	<input type="text" value="1"/>	<input type="text" value="11.072+1.456"/>	<input type="text" value="1.353+0.763"/>	<input type="text" value="1.391+0.723"/>
<input type="text" value="2015-07-01 13:31:24"/>	<input type="text" value="1"/>	<input type="text" value="49.737+1.809"/>	<input type="text" value="8.232+0.910"/>	<input type="text" value="10.329+0.963"/>
<input type="text" value="YYYY-MM-DD HH:mm:ss"/>	<input type="text" value="1.00"/>	<input type="text" value="0.00 ± 0.00"/>	<input type="text" value="0.00 ± 0.00"/>	<input type="text" value="0.00 ± 0.00"/>

Series Name:

Date-Time	Dilution	133 Xe	131m Xe	133m Xe
<input type="text" value="2015-07-01 13:45:10"/>	<input type="text" value="1"/>	<input type="text" value="4.202+0.161"/>	<input type="text" value="0.294+0.104"/>	<input type="text" value="0.180+0.097"/>
<input type="text" value="YYYY-MM-DD HH:mm:ss"/>	<input type="text" value="1.00"/>	<input type="text" value="0.00 ± 0.00"/>	<input type="text" value="0.00 ± 0.00"/>	<input type="text" value="0.00 ± 0.00"/>

[+ add series](#)

Figure 5-3. Series parameters

After entering sample data, there are a number of models that can be selected to compare against the entered data. The dropdowns allow the user to select a model. The magnitude of the model can be adjusted by use of a Scalar value for comparison with the entered data. This will be apparent on charts with an axis in the Activity Domain. The start time of the model can be adjusted to match with charts that are measured with Time.


Models

Model Type	Irradiation Time	Separation Time	Scalar	Start Date-Time		
<input type="text" value="bwr"/>	<input type="text" value="1260 days"/>	<input type="text" value="infinity"/>	<input type="text" value="1"/>	<input type="text" value="2015/06/08 00:00:00"/>	<input type="text" value="x"/>	<input type="text" value="x"/>
<input type="text" value="mipf"/>	<input type="text" value="5 days"/>	<input type="text" value="3 hours"/>	<input type="text" value="1"/>	<input type="text" value="2015/06/08 00:00:00"/>	<input type="text" value="x"/>	<input type="text" value="x"/>
<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value=""/>	<input type="text" value="1.00"/>	<input type="text" value="YYYY/MM/DD HH:mm:ss"/>	<input type="text" value="x"/>	

Figure 5-4. Model selection and parameters

After the data has been entered and models chosen, charts can be created to confirm the models represent the dataset. It is important to note that while Ratio vs Ratio plots may match a model, it is possible that Activity vs Time or Activity vs Activity may not match a model. When that occurs, it implies the model does not represent the entered data. Uncertainties will be displayed on the chart as black bars when inputted in the sample data.

Chart



[+ create chart](#)

Figure 5-5. Isotope Analyzer plotting




Figure 5-6. Multi-isotope plot

6.0 Station Health

6.1 Station Health Application

The Station Health application is an overview of station health categorized by the measurement type (Sample PHD, Gas Background PHD, and Quality Check PHD). This application displays tabular information of each station's health based on their measurements. These health metrics per measurements are cc of xenon, quality indicator (Good, Fair, or Bad status), Last Full Sample PHD received date, number of measurements received in the last seven days, measurement reception interval (time span based), indication of past due measurement (red is past due, green is on time), and the time span of the past due measurement. For details, see Figure 5.1.

Station Health

Station	PHD Type	Measurement	Last Full Sample PHD Received	Number received in last 7 days	Interval	Past Due	Overdue
ALU-001	SAMPLEPHD		2017/01/12 21:47:47 UTC	0			
	GASBKPHD		2017/01/12 21:47:47 UTC	0			
	QC PHD		2017/01/12 21:47:47 UTC	0			
ALU-002	SAMPLEPHD		2017/01/12 21:47:47 UTC	0			
	GASBKPHD		2017/01/12 21:47:47 UTC	0			
	QC PHD		2017/01/12 21:47:47 UTC	0			
ALU-003	SAMPLEPHD		2017/01/12 21:47:47 UTC	0			
	GASBKPHD		2017/01/12 21:47:47 UTC	0			
	QC PHD		2017/01/12 21:47:47 UTC	0			
ALU-004	SAMPLEPHD		2014/03/03 13:53:28 UTC	0	a day		3 years
	GASBKPHD		2014/03/02 19:46:06 UTC	0	a day		3 years
	QC PHD		2014/03/03 14:30:58 UTC	0	a day		3 years
ALU-005	SAMPLEPHD		2017/01/12 21:47:47 UTC	0			
	GASBKPHD		2017/01/12 21:47:47 UTC	0			
	QC PHD		2017/01/12 21:47:47 UTC	0			
ALU-006	SAMPLEPHD		2017/01/12 21:47:47 UTC	0			
	GASBKPHD		2017/01/12 21:47:47 UTC	0			
	QC PHD		2017/01/12 21:47:47 UTC	0			
ALU-007	SAMPLEPHD		2017/01/12 21:47:47 UTC	0			
	GASBKPHD		2017/01/12 21:47:47 UTC	0			
	QC PHD		2017/01/12 21:47:47 UTC	0			
ALU-008	SAMPLEPHD		2017/01/12 21:47:47 UTC	0			
	GASBKPHD		2017/01/12 21:47:47 UTC	0			
	QC PHD		2017/01/12 21:47:47 UTC	0			
ALU-009	SAMPLEPHD		2017/01/12 21:47:47 UTC	0			
	GASBKPHD		2017/01/12 21:47:47 UTC	0			
	QC PHD		2017/01/12 21:47:47 UTC	0			
ALU-010	SAMPLEPHD		2017/01/12 21:47:47 UTC	0			
	GASBKPHD		2017/01/12 21:47:47 UTC	0			
	QC PHD		2017/01/12 21:47:47 UTC	0			
ALU-011	SAMPLEPHD		2017/01/12 21:47:47 UTC	0			
	GASBKPHD		2017/01/12 21:47:47 UTC	0			
	QC PHD		2017/01/12 21:47:47 UTC	0			
ALU-012	SAMPLEPHD		2017/01/12 21:47:47 UTC	0			
	GASBKPHD		2017/01/12 21:47:47 UTC	0			
	QC PHD		2017/01/12 21:47:47 UTC	0			
ALU-013	SAMPLEPHD	ccXe	2017/01/12 18:49:59 UTC	14	a day		
	GASBKPHD		2017/01/12 18:48:56 UTC	14	a day		
	QC PHD		2017/01/12 19:34:50 UTC	28	12 hours		
ALU-014	SAMPLEPHD	ccXe	2017/01/12 18:44:28 UTC	14	a day		
	GASBKPHD		2017/01/12 18:44:07 UTC	14	a day		
	QC PHD		2017/01/12 18:44:07 UTC	14	a day		

Figure 6-1. Station Health Application

7.0 Live View

The Live View application (Figure 6.1) is located in the menu next to the Stations link within the RN Analyst role. Live View provides a current display of the latest 10 measurements sent to Watchmen. This screen refreshes every 10 seconds.

Each measurement has some metadata along with a spectra thumbnail. You can hover over the spectra just as in the Search applications to get extra details (Figure 2.1). The Live View also displays the status indicators (Section 3.1.7). You can click on the spectra to navigate to the measurement in the Search application (Section 3.0).




Figure 7-1. Live View Application

8.0 SOH Overview

8.1 Station Performance Overview

The Station Performance overview screen gives an overall view of Watchmen’s stations and their performance. The overview is based on a historical view of quality indicators. The three statuses for any indicator are Bad, Fair, and Good. The status of the indicators displays as the worst-performing indicator in any given time frame—that is an indicator that is in Bad state will override a Good or Fair indicator. By default, the overview screen displays all stations’ performance over the past seven days with the stations in Bad performance ordered at the top of the screen. Each station displays a time-based horizontal band per detector. A station’s current performance is based on the most recent quality indicator’s status. A user can optionally change the time frame of a station to get historical performance. A user can also re-sort the stations based on status by clicking on the status legend at the top of the screen. The status legend also shows the count in each status marker (see Figure 7.1)

Upon clicking on a station graph, the application navigates to the station detail screen.





Figure 8-1. Station Performance

8.2 Station Health

On the station detail screen (see Section 9.1 for full station detail documentation), the heading “Station Health” contains indicator information that defines this station’s health. Similar to the Station Performance Overview graphs (see Figure 7.1), the Station Health section displays a set of timebased horizontal band of quality indicators. These quality indicators are rolled up to the topmost indicator, which is the Station Overview. The Station Overview is equivalent to the Station Performance Overview station graph. Hovering over any point in an individual indicator gives a detail of the indicator value and the details of the status—Bad, Fair, Good (see Figure 7.3). Clicking on any point in an individual indicator will navigate to the containing measurement (Figure 3.6). White space indicates no data are available at a specific time because of expected or unexpected missing data. Each individual indicator’s data for the given time range can be exported to CSV via the icon to the right of the indicator name (). See Figure 7.2.




Figure 8-2. Station Health





Figure 8-3. Station Health Hover

9.0 Ad Hoc Query Tool

The Ad Hoc Query Tool allows the user to compare data from different stations and detectors. Data from the same detector but different variables can also be compared. Upon entering the Ad Hoc Query Tool page, the user will be met with three empty fields to select the Station, Detector and Variable to plot. There is also a date range selector on the right-hand side that defaults to the range of a year.



The initial Ad Hoc Query Tool page features a search interface at the top. It includes three input fields labeled 'Station', 'Detector', and 'Variable', each with a 'Search...' placeholder. To the right of these fields is an 'Add to Graph' button. Further right is a 'Date Range' selector showing a default range from '01/28/2021 19:49' to '01/28/2022 19:49', with 'x' icons for clearing the selection. Below the search fields is a table with headers 'Station', 'Detector', and 'Variable'. The main area of the page is a large rectangle with a magnifying glass icon and the text 'No data to graph'. Below this text is a small instruction: 'Select a station, detector and variable and press the 'Add to Graph' button. The first two selections will be shown on the graph.'

Figure 9-1. Initial Ad Hoc Query Tool Page

Once all the fields are filled out the user can then press the "Add to Graph" button. If the data exists it will be plotted on a graph.




Figure 9-2. Plotting a Single Variable

The user can plot data from two different variables at a time. Additional variables beyond the first two on the table will not be shown on the graph. Variables can be removed from the graph by pressing the "trash" button to remove that entry from the table.



Figure 9-3. Plotting Multiple Variables

Clicking a point on the graph will open a link to the associated measurement in a new tab if a measurement exists for that point. The date range can also be adjusted at any time to update the time span for the data shown on the graph.




Figure 9-4. Adjusting the Date Range

The user can toggle between a line or scatter series for the chart by pressing the button on the top left of the graph. They can also use the download button on the bottom left to export the data to a CSV.




Figure 9-5. Line Series Toggled

There is a moving range chart below the main chart which shows the absolute difference between each value compared to the previous value.



Figure 9-6. Moving Range Chart

10.0 World

The World page displays where stations are located around the world along with approximate locations of nuclear and medical facilities that may emit radiation. Each type of site is displayed on its own layer. Clicking on the name of a layer in the box at the top right of the screen toggles that layer's visibility. Clicking on a marker will display information about that station or facility. Note that some facilities are very close to each other, and it may be difficult to distinguish them from one another when zoomed far out.

A station's location must be updated manually by a user with administrator privileges. This can be done by navigating to the Stations page, selecting a station, clicking the Edit Station button, and adjusting the Latitude and Longitude fields.




Figure 10-1. World Page

11.0 Administration

The administrative functions previously contained in a separate interface are now included in the Navigation. To open, select the administrative option of interest from the menu. The options are:

- Pending Reviews
- Stations
- Indicator Configuration.
- Users

11.1 Pending Reviews

A station may only be assigned to one reviewer. The stations should be assigned to users with a role of Analyst. If a change in station assignments is made, any samples loaded for stations assigned to users after a change has been made will start appearing in the Pending Reviews application for the changed user. However, already assigned reviews will not be affected.

To change a station's measurement review to a different user, select the desired Station and 'Unreviewed' Quality. Click on the row corresponding to the measurement that you are reassigning and select the Reassign button. In the pop-up menu, select the new user from the drop-down menu. To highlight multiple samples, hold down the Shift or Ctrl key and click to select multiple rows (Figure 9.1).

NOTE: You will see the same filters as the Search application in the Get Measurements panel (see Section 3.1). You can use these to filter the measurements to assign.




Figure 11-1. Pending Sample Reviews Application

11.2 Stations

The Stations application contains a list of all the stations currently available in the database for viewing in the Watchmen interface. The list displays the Station Code, Country Code and Station Type, where are also the categories by which the list may be filtered.

Select Station	Country Code	Station Type	
<input type="text" value="Search..."/>	<input type="text" value="Search..."/>	<input type="text" value="Search..."/>	<input type="button" value="New Station"/>
Station	Country Code	Station Type	
ARB15	US	RASA	
ARP01	US	RASA	
ARP03	AR	RASA	
ARX01	xx	Not Set	
AUP04	AU	RASA	
AUP05	AU	RASA	
AUP06	AU	RASA	

Figure 11-2. Station List

In the upper right, there is a button to add a New Station. The new station dialog designates which attributes are required with a red star to the right of the field. Once completed, click Done will open the Station's Information panel (which will mostly be void, because no data will have been loaded to the database for it). To return to the Station list, click the X in the select Station field to remove the filter.

New Station

Information

Station Name

Country Code

Latitude

Longitude

Elevation

POCID

Station Type

Station Model

Calculation Type

Type

Status

Description

Assigned To

Indicator Configuration

Figure 11-3. New Station Dialog

Select Station

EG001

New Station

Station EG001

Information

Country CodeLatitude 0Longitude 0Elevation 0POCID 0

Station TypeBetaGammaStation ModelNot SetCalculation TypePNNL BetaGamma 7 ROITypeStatus

DescriptionTRIOAssigned To1122TeStErIndicator Configuration

Edit StationDelete StationManage Station Models

Detectors

No detector data for station EG001

Trend Range

Date Range01/08/2022 00:4702/08/2022 00:47

CalibrationNo Calibrations

Open Samples ReportDownload Samples Report

No trend data to display

Station Health

Quality Status

Figure 11-4. New Station Information Panel

Clicking on the Stations Code from the Station list will open the Information panel for that station.




Figure 11-5. Station Information Panel

Clicking on the Edit button for a detector opens the Detector Calibrations page in a separate tab.

Station USX75

Detector USX75_001

Expected Gain Values

Beta channel: Gamma channel:

Custom Calibrations

Name	Valid From	Valid To			
RL16 2016 MITER	06/03/2019 00:00	01/28/2022 00:00	<input type="button" value="Edit"/>	<input type="button" value="View"/>	<input type="button" value="Delete"/>

Original Calibrations

ID	Calibration Date	Sample Start	Sample Stop		
1862	05/15/2007 10:33	01/26/2013 17:04	02/17/2013 19:02	<input type="button" value="Clone"/>	<input type="button" value="View"/>
812	05/15/2007 03:33	06/09/2012 13:40	12/07/2014 15:53	<input type="button" value="Clone"/>	<input type="button" value="View"/>
721	05/15/2007 03:33	08/04/2011 01:43	06/09/2012 02:10	<input type="button" value="Clone"/>	<input type="button" value="View"/>
8	05/15/2007 10:33	05/10/2007 13:10	06/14/2011 20:02	<input type="button" value="Clone"/>	<input type="button" value="View"/>
841	06/01/2012 03:33	05/10/2007 06:11	05/10/2007 06:11	<input type="button" value="Clone"/>	<input type="button" value="View"/>
281	05/24/2006 10:38	07/26/2006 08:22	07/26/2006 08:22	<input type="button" value="Clone"/>	<input type="button" value="View"/>

Figure 11-6. Detector Calibrations

On the Detector Calibrations page, expected gain values can be adjusted for the detector, and the detector's calibrations can be viewed, added, and cloned.

Clicking the Add Calibration button allows a user to add a new custom calibration from the calibration of a selected sample file. The name, valid from, and valid to fields for a custom calibration can be changed by clicking on the Edit button for the calibration.

Users with administrative permissions can remove custom calibrations for a detector by clicking the Delete button for that calibration. A calibration's data can be viewed by clicking the View button for the calibration.




Figure 11-7. Viewing Calibrations

A scrollable window opens, allowing the user to view the ROI Limits, Ratios, Gamma, Beta, Gamma Beta Efficiency, and Software tables for a calibration.

The Information panel contains the following information:

- Station Code – A code given to each station that is the list of stations displayed throughout the interface
- Country Code – The two-letter code of the country operating the station
- Station Type – A feature not currently implemented in this software
 - Not Set (no calculations will be performed)
 - PNNL BetaGamma 7 ROI
 - FOI BetaGamma 10 ROI
 - CEA Basic SPALAX
 - Melusine

To manage the station information, click the Edit Station button at the bottom of the Information panel. This will enable in-line editing. Make the desired changes, and then select the Done checkmark button to save the changes or the X Cancel button to cancel.

The Detectors information provides energy calibration information. The information are coefficients of a polynomial fit function of the form Ax^2+Bx+C where A , B , and C are numbers provided in the chart on the page. There are two type of calibration information provided for each the beta detector and the gamma detector: Channel to energy (C2E) coefficients and energy to channel (E2C). These coefficients are used when performing analysis. They can be updated manually by clicking the edit icon to the right of the coefficients. Measurements will not be reanalyzed automatically when these coefficients are updated.

For station trend plotting, select the Station name from the drop-down menu. If trend data exists for the station, it is automatically plotted in the ‘Trend Range’ panel. The default date range is one month. The date range can be modified, and the chart will automatically update. If the station has custom calibrations, they will appear in the ‘Calibration’ dropdown. The trend chart will display the data for the selected calibration.




Figure 11-8. Station Trend Range Panel

Clicking the ‘Open Samples Report’ button in the Trend Range will open a pop-up with tabular data for the selected calibration. You can also click the ‘Download Samples Report’ button to get a Comma Separated Values (CSV) export of the trend data.

NOTE: The “Station Health” section is detailed in Section 8.2.

11.3 Detector Configuration

Clicking a detector link from the stations page will open the detector configuration dialog.



Station **USX75**
Detector **USX75_008**
Custom Calibrations
[Add Calibration](#)

Figure 11-9. Default Detector Configuration

From this dialog a user may upload alternate detector calibrations (from that which is included in the measurement files themselves) for performing additional analysis. Clicking the “Add Calibration” allows the user to select a measurement file containing the custom calibration, set a descriptive name, and the period for which the custom calibration is valid (“Valid Start” is required, “Valid Stop” is optional). The default for “Valid Start” is the current date/time, but may be changed.




Station **USX75**
Detector **USX75_008**
Custom Calibrations

Sample File	Name	Valid start	Valid stop
Choose sample file... Browse		09/21/2021 09:14:02	

☒

Figure 11-10. Adding Custom Calibration

If a calibration from a different system and/or detector is attempting to be used, a dialog will ask the user to verify that they still want to import the mis-matched calibration. Selecting “Ok” will close the dialog.



Station Detector Warning

The file contains a different station or detector than the station and detector you're uploading the calibration to.

Station	File Station
USX75	USL16
Detector	File Detector
USX75_008	USL16_005

Click **OK** to use the calibration in this file for station **USX75** and detector **USX75_008**.

Figure 11-11. Mis-match Detector Dialog

Once required entries are filled, the checkmark button is enabled. After clicking the checkmark button, the file is uploaded, parsed, inserted into the database, and then page is updated.

Station **USX75**
 Detector **USX75_008**
 Custom Calibrations

[Add Calibration](#)

Name	Valid Start	Valid Stop			
Alt Calibration 2	06/01/2018 10:23:59	09/08/2021 17:00:00	Edit	View	Delete
Alt Calibration	09/13/2018 14:30:36		Edit	View	Delete

Figure 11-12. **Uploaded Detector Calibrations**

After uploading a custom calibration, when measurement files are subsequently analyzed (either upon imported, or user triggered recalculate), if the measurement’s acquisition start time falls within the custom calibration’s valid period, then additional analysis will be performed on the measurement with the custom calibration and its results may be viewed on the Sample Review Panel but selecting the calibration from the Calibration dropdown select.

11.4 Indicator Configuration

Users can administer and change the configuration of limits for a station’s quality indicators (see Section 7.2 for a definition of indicators). Start by going to the “Indicator Config” menu item (found in the Administrator tab on the header). You can then modify a configuration by selecting the configuration from the “Select a Config” drop-down menu. Or you can start a new configuration from a pre-defined system template by clicking “Start from Template” (see Figure 9.3). Selecting an existing config navigates to the config detail screen with the selected configuration populated. Selecting “Start from Template” navigates to the same config detail screen with a default set of values populated.




Figure 11-13. **Indicator Administration Start**

On the Station Indicator Administration detail page users can create or modify station configurations. Station configurations can be assigned to stations on the Station Detail screen (see Section 9.1).

The Station Indicator Administration detail page can set the configuration name, along with details for each quality indicator. For each quality indicator, users can adjust the limit label, and adjust the limit value. If the quality indicator value crosses a limit value, the indicator’s status is assigned within that limit. See Figure 9.4 for an overview of the screen.

Users can also toggle email notifications for an indicator’s limits. If a quality indicator value crosses a limit and the email notification is enabled for that limit, a notification is sent to the configured email address at the time of calculation. See Figure 9.5.

NOTE: A weekly summary of all indicators is sent regardless of notification settings in the configurations.

Once satisfied with the configuration, users can save the new configuration—or overwrite an existing configuration—by clicking the appropriate button (located under the “Config Name” field). For a description of an indicator configuration, see Figure 9.6.

[← Return to Configuration Selection](#)

Station Indicator Administration

Last Updated: N/A

Config Name:

[Save As New Configuration](#)

Acquisition

BAD	<input type="checkbox"/>	<input type="text"/>	<input type="text" value="41004"/>	\$
WARN	<input type="checkbox"/>	1% to 2%	<input type="text" value="40602"/>	\$
GOOD	<input type="checkbox"/>	<input type="text"/>	<input type="text" value="39798"/>	\$
WARN	<input type="checkbox"/>	-1% to -2%	<input type="text" value="39396"/>	\$
BAD	<input type="checkbox"/>	<input type="text"/>		

Air

GOOD	<input type="checkbox"/>	<input type="text" value="≥ 0"/>	<input type="text" value="0"/>	m³
BAD	<input type="checkbox"/>	<input type="text" value="< 0"/>		

ccXe

BAD	<input type="checkbox"/>	<input type="text" value="> 1.5"/>	<input type="text" value="1.5"/>	cc
WARN	<input type="checkbox"/>	<input type="text"/>	<input type="text" value="1.25"/>	cc
GOOD	<input type="checkbox"/>	<input type="text"/>	<input type="text" value="0.75"/>	cc
WARN	<input type="checkbox"/>	<input type="text"/>	<input type="text" value="0.5"/>	cc
BAD	<input type="checkbox"/>	<input type="text" value="< 0.5"/>		

Collection

BAD	<input type="checkbox"/>	<input type="text"/>	<input type="text" value="44064"/>	\$
WARN	<input type="checkbox"/>	1% to 2%		

Figure 11-14. Station Indicator Administration Configuration

Acquisition

BAD	<input type="checkbox"/>		
WARN	<input checked="" type="checkbox"/>		41004 \$
GOOD	<input type="checkbox"/>		40602 \$
WARN	<input type="checkbox"/>	-1% to -2%	39798 \$
BAD	<input type="checkbox"/>		39396 \$

Figure 11-15. Toggle E-mail Notifications



Acquisition			
Indicator Status	BAD	<input checked="" type="checkbox"/>	
Toggle E-Mail	WARN	<input type="checkbox"/>	1% to 2%
Limit Label	GOOD	<input type="checkbox"/>	
Limit Value	WARN	<input type="checkbox"/>	-1% to -2%
	BAD	<input checked="" type="checkbox"/>	

Figure 11-16. Indicator Configuration Detail

11.5 Users

The Manage Users screen is used to assign roles to users which controls which review queues they have access to.

Manage Users

<div>New User</div>							
User		Enabled	Analyst	Lead Analyst	Snr. Analyst	Process Eng.	Dir. of Ops.
Admin		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
DisabledUser		<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
NewUser		<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>


NOTE: The new user button will be removed in a future version.

The first column shows common name of the user, hovering the mouse over the name shows the full name of the user if it is different.

The checkboxes can be used to select the roles (Analyst, Lead Analyst, Senior Analyst, Processing Engineer, or Director of Operations) a user is assigned. The Analyst role will be automatically added to a user granted the role(s) Lead Analyst, Senior Analyst, Processing Engineer, or Director of Operations. The Lead Analyst role will be automatically added to a user granted the Director of Operations role. Changes are saved automatically. Inherited roles will not be unchecked along with the heir role (e.g. unchecking Lead Analyst role will not uncheck the Analyst role). The “Enabled” role controls access to Watchmen; any user not enabled will be shown an error when trying to access Watchmen. A user is not assigned a role, besides “Enabled,” will be able to view measurements, but will not be capable of reviewing them. To remove an existing user, click the trash can to the right of the username. Only users without references (no assigned or reviewed measurements) in the Watchmen application can be deleted.

New users added via New User button, or through the web request, will not be enabled, nor will they be assigned roles automatically. These new users will be filtered to the top of the users list. Additionally, a bell icon will be added to the Administrator’s menu, and the Users menu item to notify Administrative users that a user is awaiting a role. Enabling the user, or assigning them roles, will turn off the notification, however while on the Users page, the new user will remain at the top, as there may be several roles that are needed to be applied to the new user.

12.0 Information

By clicking on the question mark in the menu bar () , users can access the Information screen. This screen displays information about the Watchmen application. The information sections are defined in the following areas. See Figure 10.1 for an illustration of the information screen.

Server Information

Components	Environment
War	Server Apache Tomcat/7.0.35
Version 2.3	Java VM Name IBM J9 VM
Built 2017/01/24 20:36:40	Java VM Vendor IBM Corporation
	Java VM Version 1.7.0
	OS Name Linux
	OS Version 2.6.32-131.0.15.el6.x86_64
	Max Memory 1024 MiB
	Total Memory 1024 MiB
	Free Memory 317.4 MiB
	Processors 8

Calculation Information

[10 Region Calculation Definition Report](#)

Legal

This computer software was prepared by Battelle Memorial Institute, hereinafter the Contractor, under Contract No. DE-AC05-76RL0 1830 with the Department of Energy (DOE). All rights in the computer software are reserved by DOE on behalf of the United States Government and the Contractor as provided in the Contract. You are authorized to use this computer software for Governmental purposes but it is not to be released or distributed to the public.

This material was prepared as an account of work sponsored by an agency of the United States Government. Neither the United States Government nor the United States Department of Energy, nor the Contractor, nor any of their employees, nor any jurisdiction or organization that has cooperated in the development of these materials, **makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness or any information, apparatus, product, software, or process disclosed, or represents that its use would not infringe privately owned rights.**

Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or any agency thereof, or Battelle Memorial Institute. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or any agency thereof.

PACIFIC NORTHWEST NATIONAL LABORATORY
operated by
BATTELLE
for the
UNITED STATES DEPARTMENT OF ENERGY
under Contract DE-AC05-76RL01830

Notice

This software bundles the following NPM packages:

backbone	https://spdx.org/licenses/MIT
backbone-react-component	MIT (https://github.com/magalhas/backbone-react-component/blob/master/LICENSE-MIT)
bootstrap	MIT (https://github.com/twbs/bootstrap/blob/master/LICENSE)
d3	BSD (https://github.com/mbostock/d3/blob/master/LICENSE)
filesaver.js	https://spdx.org/licenses/MIT
flux	https://spdx.org/licenses/BSD-2-Clause
font-awesome	OFL-1.1 (http://scripts.sil.org/OFL), MIT (http://opensource.org/licenses/mit-lic)
jquery	MIT (https://github.com/jquery/jquery/blob/2.1.4/MIT-LICENSE.txt)
keymirror	Apache-2.0 (http://www.apache.org/licenses/LICENSE-2.0)
leaflet	https://spdx.org/licenses/BSD-2-Clause

Figure 12-1. Information Detail

12.1 Server Information

This identifies information about the web server hardware and versioning, as well as the software build date.

12.2 Calculation Information

This is a link to a PDF document that illustrates the 10-region calculation process for Watchmen.

12.3 Legal

This section is a legal software notice from PNNL operated by Battelle, sponsored by the U.S. Department of Energy.

12.4 Notice

This section reveals all third-party software licenses used in the Watchmen application.

12.5 Developer Information

This section contains a link to the development REST API for documentation.

Pacific Northwest National Laboratory

902 Battelle Boulevard
P.O. Box 999
Richland, WA 99354
1-888-375-PNNL (7665)

www.pnnl.gov