PNNL-23832



Toxicity of Uranium Adsorbent Materials using the Microtox Toxicity Test

DOE-NE Milestone Report M3FT-15PN0310056

Jiyeon Park Robert Jeters Gary Gill Li-Jung Kuo George Bonheyo

October 2014



Proudly Operated by Battelle Since 1965

DISCLAIMER

This report was prepared as an account of work sponsored by an agency of the United States Government. Neither the United States Government nor any agency thereof, nor Battelle Memorial Institute, nor any of their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or any agency thereof, or Battelle Memorial Institute. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or any agency thereof.

PACIFIC NORTHWEST NATIONAL LABORATORY

operated by

BATTELLE


for the

UNITED STATES DEPARTMENT OF ENERGY under Contract DE-AC05-76RL01830

Printed in the United States of America

Available to DOE and DOE contractors from the Office of Scientific and Technical Information, P.O. Box 62, Oak Ridge, TN 37831-0062; ph: (865) 576-8401 fax: (865) 576-5728 email: reports@adonis.osti.gov

Available to the public from the National Technical Information Service 5301 Shawnee Rd., Alexandria, VA 22312 ph: (800) 553-NTIS (6847) email: order@ntis.gov order@ntis.gov order@ntis.gov order@ntis.gov order@ntis.gov <a href="mailto:http://www.ntis.g



Toxicity of Uranium Adsorbent Materials using the Microtox Toxicity Test

DOE-NE Milestone Report M3FT-15PN0310056

Jiyeon Park, Robert Jeters, Gary Gill, Li-Jung Kuo, and George Bonheyo

Marine Sciences Laboratory Pacific Northwest National Laboratory 1529 W. Sequim Bay Road Sequim, WA 98382

October 2014

Prepared for the U.S. Department of Energy under Contract DE-AC05-76RL01830

Pacific Northwest National Laboratory Richland, Washington 99352

Executive Summary

The Marine Sciences Laboratory at the Pacific Northwest National Laboratory evaluated the toxicity of a diverse range of natural and synthetic materials used to extract uranium from seawater. The uranium adsorbent materials are being developed as part of the U. S. Department of Energy, Office of Nuclear Energy, Fuel Resources Program. The goal of this effort was to identify whether deployment of a farm of these materials into the marine environment would have any toxic effects on marine organisms. Toxicity was evaluated using the well-recognized and accepted Microtox® toxicity test. The Microtox® assay uses a microorganism found in seawater that is susceptible to a wide range of toxins. This established toxicity assay is quick, reliable and cost-efficient.

Direct contact with 5 out of 10 adsorbents were toxic (reached EC50), but only when the adsorbent material was at very high concentrations. No toxicity was observed with seawater that was in contact with the absorbent material from 63 different adsorbent materials. Based on testing to date, there does not appear to be any significant concerns with toxicity of adsorbent materials; Extraction from seawater could be performed with minimal impact to marine environment.

Materials and Methods

The toxicity of adsorbent materials were evaluated using the Microtox® assay to determine if the supplied adsorbent materials and seawater effluent from uranium extraction columns have any toxic effects on marine organisms, The assay measures the luminescence of the bioluminescent marine bacterium *Aliivibrio fischeri* (ATCC 49387), an indicator of metabolism in the organism, after a 30-minute (1) exposure to the sample material. The Microtox® assay is a useful tool to evaluate the toxicity of a wide range of metals (2, 3) and organic (4, 5) compounds that are either in liquid or solid phase. In addition, the Microtox® test is low-cost and easy to operate and has a long history of use for determining toxicity in environmental samples (6-8). The Advanced Monitoring Systems (AMS) Center that is part of the Environmental Technology Verification (ETV) Program has evaluated, tested and validated the Microtox® assay for assessing various contaminants in water samples. The ETV is supported by the U.S. Environmental Protection Agency (EPA) (9). The Microtox® assay has also been used by U.S. Geological Survey and National Oceanic and Atmospheric Administration (NOAA) to assay sediment from the Delaware Bay and surrounding areas for toxicity (10).

Aliivibrio fischeri culture was grown overnight in ALNa Broth (1) at 22°C with vigorous shaking. The overnight culture was used to inoculate 10 ml of ALNa broth, and the culture was incubated with vigorous shaking at 22°C . Cells were grown until the optical density at 590 nm (OD_{590}) was 0.25. Cells were harvested by centrifugation at 4150 rpm for 5 min followed by washing in 10ml 3% NaCl; pH 7. Cells were then resuspended in 10 mL 3% NaCl at pH 7. The cell suspension was mixed with test materials for 30 minutes, and the luminescence was measured using a Synergy HT microplate reader (Biotek, Winooski, VT, USA).

The decrease in bacterial luminescence (INH%) (11) was calculated as shown below. If the luminescence decreased more than 50% after exposure to sample materials for 30 minutes, the concentration of toxicant is designated as the effective concentration (EC $_{50}$) (11). Values obtained with cells only control sample (suspended in 3% NaCl or seawater) was set to 100%, and the luminescence of each sample was represented as a percentage of cells only value. ZnSO $_{4}$ (12) was used as control toxicant (0.625, 1.25, 2.5, and 5 μ g/mL to validate the Microtox® assay results.

The decrease in bacterial luminescence (INH%)

INH% =
$$100 - \left(\frac{ITT}{IT0 \times KF}\right) \times 100$$

$$KF = \frac{ICT}{IC0}$$

IC0: luminescence of control at t = 0

ICT: luminescence of control after 'T' minutes exposure time

IT0: luminescence of test sample at t = 0

ITT: luminescence of test sample after 'T' minutes exposure time

KF: correction factor based on the control/blank

A positive response means a decrease in luminescence, while a negative response means the luminescence increased.

EC50 (50% Effective Concentration)


The concentration of the toxicant (mg/L) that causes a 50% reduction in light after exposure for 'T' minutes

Results and Discussion

1. Toxicity of adsorbent materials via direct contact

Supplied adsorbent materials were tested for direct contact toxicity using increased concentrations. Due to limited number of adsorbent materials currently available, only 10 samples were tested, and hence the results for direct testing should be considered preliminary.


- 1) Oak Ridge National Laboratory (Richard Mayes)
 - YF130313C KOH Condition did not have an EC₅₀ up to 25% v/v.
 - YF130226C Copolymer-PAN had an EC₅₀ at 25mg/ml.




• The Microtox® assays on YF130313C KOH Condition and YF130226C Copolymer-PAN were done one time due to the limited amount of supplied sample materials.

2) University of North Carolina (Carter Abney/Lin group)

- CA-4139: Porous Titanium Dioxide had an EC₅₀ at 1.1mg/ml.
- CA-4149: Porous Titanium Phosphate did not have an EC₅₀ up to 6.4mg/ml.
- MC-Precursor Unfunctionalized Mesoporous Carbon did not have an EC₅₀ up to 10mg/ml.
- MC-Phosphate (612) Phosphoryl-functionalized Mesoporous Carbon had an EC₅₀ at 5mg/ml.




• CA-4139: Porous Titanium Dioxide and CA-4149: Porous Titanium Phosphate were tested in duplicated experiments. Graphs show the mean value and standard deviation is shown as error bars.




The Microtox® assays on MC-Precursor Unfunctionalized Mesoporous Carbon and MC-Phosphate Phosphoryl-functionalized Mesoporous Carbon were done one time due to the limited amount of supplied sample materials.

3) Chinese Academy of Science

- UHMWPE-g-PAPA-co-PAO #1 had an EC₅₀ at 0.625mg/ml.
- UHMWPE-g-PAPA-co-PAO #2 had an EC₅₀ up to 0.3125mg/ml.
- UHMWPE-g-PAPA-co-PAO #3 did not have an EC₅₀ up to 2.5mg/ml.
- UHMWPE-g-PAPA-co-PAO #4 did not have an EC₅₀ up to 2.5mg/ml. All fiber types 1~4 were tested in triplicate. Graphs show the mean value and standard deviation is shown as error bars.






• UHMWPE-g-PAPA-co-PAO Chinese Functionalized Fiber Type 3 and 4 showed no signs of toxicity. Increase in luminescence is likely due to compound present on functionalized fiber. Some select compounds have been shown to increase luminescence at high concentrations.

2. Toxicity of seawater effluent collected from uranium extraction columns

Toxicity was evaluated in seawater that passed through columns containing adsorbent materials (effluent test). The Microtox assay was performed with increasing amounts of seawater effluent up to 75% v/v because a minimum amount of *Aliivibrio fischeri* cells were required to measure luminescence. None of seawater effluent samples had an EC_{50} . The experiments were done in triplicate. Mean values are shown and standard deviations are indicated as error bars.


1) Oak Ridge National Laboratory

38H21: 3 daysTSU 20A: 1 weekCostas Fiber: 8 weeks




2) Japan

• Japan 2010: 3 days




3) University of Alabama

AO Fiber: 6 weeksDA Fiber: 6 weeksSS Fiber: 6 weeks



4) Chinese Academy of Science

• Chinese Fiber #2: 3 days





3. Time-course experiments to determine toxicity of seawater effluent collected from uranium extraction columns


In order to determine whether a longer exposure of the adsorbent to seawater could have toxic effects on marine organisms, time-course experiments were done on seawater effluents samples that were in contact with 55 different adsorbent materials. Seawater effluent samples were collected at three different time points; 1 hour, 1 day, and 1 week after exposure to adsorbent materials in uranium uptake columns. Samples up to 75% v/v were tested in order to measure luminescence from *Allivibrio fischeri* to


evaluate toxicity of the sample. Three independent experiments were done on these samples, and all samples up to 75% v/v did not have an EC₅₀. The experiments were done in triplicate Mean values are shown and standard deviations are indicated as error bars.


1) Oak Ridge National Laboratory (42 adsorbents)














2) The City University of New York (2 adsorbents)



3) University of Maryland (7 adsorbents)



4) Chinese Academy of Science (4 adsorbents)



Summary

We tested the toxicity of direct contact with 10 adsorbent materials at different concentrations. Due to the limited amount of supplied adsorbent samples, our results are preliminary.

Adsorbent	EC ₅₀ at	Source	
YF130313C KOH Condition	No up to 25% v/v	Oak Ridge National Laboratory	
YF130226C Copolymer-PAN	$Yes \\ \geq 25 mg/ml$	Richard Mayes	
CA-4139 Porous Titanium Dioxide	Yes ≥ 1.1mg/ml		
CA-4149 Porous Titanium Phosphate	No up to 6.4mg/ml		
MC-Precursor Unfunctionalized Mesoporous Carbon	No up to 10mg/ml	University of North Carolina Carter Abney/Lin group	
MC-Phosphate (612) Phosphoryl-functionalized Mesoporous Carbon	Yes ≥ 2.5mg/ml		
UHMWPE-g-PAPA-co-PAO Chinese functionalized fiber type 1	Yes ≥ 0.625mg/ml		
UHMWPE-g-PAPA-co-PAO Chinese functionalized fiber type 2	$Yes \ge 0.3125 mg/ml$	Chinasa Agadamy of Spignag	
UHMWPE-g-PAPA-co-PAO Chinese functionalized fiber type 3	No up to 2.5mg/ml	Chinese Academy of Science	
UHMWPE-g-PAPA-co-PAO Chinese functionalized fiber type 4	No up to 2.5mg/ml		

We also tested seawater effluent that came in contact with adsorbent materials in uranium extraction columns. The $Microtox^{@}$ assay was done at multiple concentrations and time points. All 63 samples listed below showed no toxic effects.

Adsorbent	Time Points	Source
38H21	3 Days	
TSU 20A	1 Week	
Kostas Fiber	8 Weeks	
JG159a		
JG159b		
JG164a		
JG192		
KMS-2		
TSU-27		
TSU-39C		
TSU-39D		
TSU-39E		
TSU-45H40		
38H26		
AFL		
11-651H75C		
11-105		
11-117		Oak Ridge National Laboratory
11-119	1 Hr, 1 Day, & 1 Week	
YF140115A-NH20H		
AI8 L1 R1		
AI8 L1 R2		
AI8 L2 R1		
AI8B16-TOW		
AI8B16-5PPI		
2T8 RMCJ		
AF3 RMCJ		
AF8 RMCJ		
AF160-2 RMCJ		
AF1FR2 RMCJ		
AF1FR3 RMCJ		
AF1B16-5PPI		
AF1B16-TOW		
AF1B17-5PPI		
AF1B17-25PPI		

Adsorbent	Time Points	Source	
AF1 F R2			
AF1 F R3			
AF1 L1 R1		Oak Ridge National Laboratory	
AF1 L1 R2	1 Hr, 1 Day, & 1 Week		
AF1 L1 R3			
AF1 L2 R1			
AF1 L2 R1 (Fiber)			
AF1 L2 R2 (Fiber)			
AF1 L2 R3 (Fiber)			
AF3			
AO Fiber			
DA Fiber	6 Weeks	University of Alabama	
SS Fiber			
RJ-02-76-CUNY	1 Hr, 1 Day, & 1 Week	The City University of	
RJ-03-66-CUNY	1 m, 1 Day, & 1 week	New York	
B2MP-G-169			
B2MP-G-180		University of Maryland	
B2MP-G-182	1 Hr, 1 Day, & 1 Week		
B2MP-G-184			
B2MP-G-185			
B2MP-G-186			
Recycled B2MP-CT-126			
Japan 2010	3 Days	Japan Atomic Energy Agency	
Chinese Fiber #2	3 Days		
UHMWPE-g-PAPA-co-PAO 1			
UHMWPE-g-PAPA-co-PAO 2	1 Hr, 1 Day, & 1 Week	Chinese Academy of Science	
UHMWPE-g-PAPA-co-PAO 3			
UHMWPE-g-PAPA-co-PAO 4			

Conclusion

Uranium adsorbent materials had little or no toxicity even when assayed with direct contact at pp-thousand levels. It is important to note that these adsorbent materials will be exposed in the marine environment at pp-billion or lower levels. Hence these results provide preliminary evidence that uranium extraction from seawater could be performed with minimal impact to marine fauna. We conclude that the adsorbent materials used to extract uranium from seawater are relatively non-toxic and likely to have minimal contact with marine organisms on a large scale

References

- 1. Girotti SL, Bolelli L, Roda A, Gentilomi G, & Musiani M (2002) Improved detection of toxic checmicals using bioluminescent bacteria. *Analytical Chimica Acta* 473:113-120.
- Fulladosa E, Murat JC, & Villaescusa I (2005) Study on the toxicity of binary equitoxic mixtures of metals using the luminescent bacteria Vibrio fischeri as a biological target. *Chemosphere* 58(5):551-557.
- 3. Tsiridis V, *et al.* (2006) Interactive toxic effects of heavy metals and humic acids on Vibrio fischeri. *Ecotoxicology and environmental safety* 63(1):158-167.
- 4. Backhaus T, Scholze M, & Grimme LH (2000) The single substance and mixture toxicity of quinolones to the bioluminescent bacterium Vibrio fischeri. *Aquat Toxicol* 49(1-2):49-61.
- 5. Lei L & Aoyama I (2010) Effect-directed investigation and interactive effect of organic toxicants in landfill leachates combining Microtox test with RP-HPLC fractionation and GC/MS analysis. *Ecotoxicology* 19(7):1268-1276.
- 6. De Zwart D & Sloof W (1983) The Microtox as an alternative assay in the acute toxicity assessment of water pollutants. *Aquatic Toxicology* 4:129-138.
- 7. Chang JC, Taylor PB, & Leach FR (1981) Use of the Microtox assay system for environmental samples. *Bulletin of Environmental Contamination and Toxicology* 26:150-156.
- 8. Ankley GT, Hoke RA, Giesy JP, & Winger PV (1989) Evaluation of the toxicity of marine sediments and dredge spoils with the Microtox bioassay. *Chemosphere* 18:2069-2075.
- 9. Anonymous (2003) ETV Joint Verification Statement. (U.S. Environmental Protection Agency and Battelle).
- 10. Station USGSBRDMER (1998) Toxicity Testing of Sediment from delaware Bay and surrounding Areas. *National Oceanic and Atmospheric Administration Office of Ocean Resources Conservation Assessment*.
- 11. Parvez S, Venkataraman C, & Mukherji S (2006) A review on advantages of implementing luminescence inhibition test (Vibrio fischeri) for acute toxicity prediction of chemicals. *Environment international* 32(2):265-268.
- 12. Dutka BJ & Kwan KK (1982) Application of four bacterial screening procedures to assess changes in the toxicity of chemicals in mixtures. *Environmental Pollution* 29:125-134.





Proudly Operated by Battelle Since 1965

902 Battelle Boulevard P.O. Box 999 Richland, WA 99352 1-888-375-PNNL (7665)

www.pnnl.gov