


Smartphone Microscope



BLOOD CELLS



SALT CRYSTALS



GRAPE LEAF



POLLEN ON FLOWER

Our 100x smartphone microscope is compatible with most phones and tablets. Experiment with the microscope and have fun observing different objects.

DIRECTIONS FOR USE



- Remove case, if needed.
- Open your camera app and choose “selfie” mode.
- Slip the Smartphone Microscope clip over your device, lining up the small glass lens with the front-facing camera.



- Slowly move the clip around until a uniform, white image appears on the screen.



- Place your camera/tablet on a flat surface in a brightly lit area.
- Place the object to view (on a slide, in a plastic baggie, or clear box) directly on the lens of the Smartphone Microscope.